

ALS-3 AUTOCLA VEVERTICAL TYPE- FULLY AUTOMATIC

HSNCODE -84192010-GST@18%

Double walled construction, Inner chamber is made up of thick stainless steel (304) & the Outer body is made up of thick Stainless Steel 304. The chamber is absolutely leak proof & can be operated at any selected point in between 5 to 20 pounds per sq. inch - with PID Digital Controller & Timer, Working Temp 121 °C to 125 °C along with S.S. basket Pressure gauge, Release cock, Safety valves (2), Drain valve (bottom), Pressure cut off switch, Solenoid Valve, Water level cut off & with foot lifting arrangement.

	Size			Rating Watts
Sr No	CM	INCH	LITRES	
A	25cmx45cm	10"x18"	22ltrs	1500
B	30cmx50cm	12"x20"	35ltrs	2000
C	35cmx55cm	14"x22"	52ltrs	3000
D	40cmx60cm	16"x24"	90ltrs	4000
E	45cmx60cm	18"x24"	95ltrs	4000
F	45cmx75cm	18"x30"	120ltrs	4500
G	55cmx75cm	22"x30"	175ltrs	5000
H	55cmx90cm	22"x36"	200ltrs	6000

OPTIONAL ACCESSORIES:

Printer interface facility with single channel logger
21 CFR software with printer interface



Image will vary; its only information