

## ALS-8 B.O.D. INCUBATOR-GMP MODEL

HSNCODE-84198960-GST@18%

Castor wheel mounted cabinet is double walled with inner chamber of stainless steel 304 & outer wall S.S. 304 with mirror finish. Double walled back of the unit is fitted with air circulation fan motor with blower for maintaining temperature uniformity throughout the chamber. Motor of 1/12HP 1440RPM. Full view inner glass door permits inspection of specimens without disturbing the same. The gap between two walls is filled with puff insulation & cooling is provided by 'Kirloskar' compressor with CFC free Unit is fitted with door operated illumination lamp inside the chamber. The front panel is provided with separate indicator lamps for main, heating & incoming voltage. Supplied complete with adjustable shelves, cord & plug to work on 220/230 volts A.C. supply. The temperature range is from 5°C to 60°C accuracy  $\pm 0.5^{\circ}\text{C}$  is controlled by Microprocessor based Digital Temperature indicator - cum - controller.

Capacity	Inner Dimension
1.5Cu.Ft(40Ltrs)	35x35x35Cm
3Cu.Ft(90Ltrs)	45x45x50Cm
6Cu.Ft(170 Ltrs)	50x50x70Cm
10Cu.Ft(280Ltrs)	55x55x90Cm
12Cu.Ft(330Ltrs)	60x60x90Cm

**Note: For BOD Incubator Stabilizer is necessary OPTIONAL**

### OPTIONAL ACCESSORIES:

RS485 Interface with software for PC Connection
Data logging with Pen-drive arrangement Single Channel
Data logging with Pen-drive arrangement 4-Channel

IMAGES:



Image will vary; its only information

SOLUTIONS (ALS)



Image will vary; its only information

AP.



*Image will vary; its only information*

ARISTOTLE LAB SOLUTIONS (ALS)