

## Thermostatic Water Bath:

This apparatus is designed for use in distilling, drying, concentrating and heat soaking chemicals or biological products.

- \* Temperature control by means of PID controller.
- \* Inner chamber with SS 304

### Technical Specification (12 Position)

|                           |                             |
|---------------------------|-----------------------------|
| Power                     | : 2000 W                    |
| No. of Opening            | : 12                        |
| Temperature-control range | : RT-99° C                  |
| Temperature-fluctuation   | : $\pm 0.2^{\circ}\text{C}$ |
| Dimension                 | : 120 (D) x 90 (H)          |
| Working Room (cm)         | : 92x32x12                  |

### Models Available:

| Model No.    | V   | Hz. | Positions |
|--------------|-----|-----|-----------|
| 1800.DNEU.01 | 230 | 50  | 6         |
| 1800.DNEU.02 | 230 | 50  | 12        |



\*\*\*